



Working Together...To Help Others





BioSkin™

Be
Better™

Bio Skin® is here to help your patients

Be Better™



www.dunbarmedical.com

Most of our products have **VIDEO APPLICATION** instructions to help you and your patients learn how to put on their brace.

See each product page at:
dunbarmedical.com and click on the "Instructions" tab to watch the video.



Table of Contents

Bio Skin Philosophy.	1
Knee	3
ROM Hinged Knees.	3
Hinged Knees.	5
Conforma™ Hinged Knees.	7
Knee Skins	9
Spine.	11
Foot & Ankle	15
Wrist	19
Elbow	21
Calf, Thigh & Hip	23
Materials.	25
Sizing Chart.Back Cover

MOTION

"A body perseveres in its state of rest, or of uniform motion, unless it is compelled to change that state by forces impressed thereon."

Sir Isaac Newton

Mathematical Principles of Natural Philosophy

The human body is designed for motion. Bones, tendons, muscles, and nerves work together to turn our bodies into complex machines capable of simple movement and remarkable feats, but sometimes injury makes a body more likely to "persevere in its state of rest."

At Bio Skin, we restore motion.

We decrease pain and enhance human performance.

We help people Be Better.







Pull-on

Front-closure



Adjustable Hinge

Gladiator™ Knee Ligament Protection

The Gladiator uses a rigid hinge to protect knee ligaments, including after reconstruction. The rigid strapping system firmly anchors the Gladiator to the leg, avoiding the need for larger, more expensive braces. Pins in the hinges may easily be set to limit flexion and/or extension.

FEATURES

- Length: 15"
- Silicone gel patella ring for compression
- Easily adjustable flexion / extension settings
- Open patella
- Pull-on **or** front-closure design
- Fits left or right

INDICATIONS

- Anterior Cruciate Ligament (ACL) Tear or Sprain
- Posterior Cruciate Ligament (PCL) Tear or Sprain
- Lateral Collateral Ligament (LCL) Tear or Sprain
- Medial Collateral Ligament (MCL) Tear or Sprain
- Knee Instability

Product

Product	Product #
Gladiator - Open Patella - Ultima™ Material	4134*
Gladiator - Open Patella - Stratus™ Material	6136*
Gladiator - Front Closure - Open Patella	5136*

Note: To indicate size, replace * with corresponding number from the sizing chart on the back cover.

Front-Closure Knee Braces

Bio Skin's front-closure knee braces use our proprietary 2SL material to create a truly customizable wrap-around design. The material itself is Velcro-compatible so the Velcro closures can adhere anywhere on the front of the brace to fit virtually any shape of leg.



Q Lok ROM™ Patella Traction System Featuring the new SlimLok™ Adjustable Hinge

The Q Lok ROM is designed to diminish the effects and aid in the rehabilitation of patellofemoral pain. The Dynamic Patella Traction System exerts a medial pull on the patella to improve tracking in the trochlear notch and increase surface contact area in the patellofemoral joint.

FEATURES

- Length: 13"
- Addresses glide and tilt
- Fits left or right
- Easily adjustable flexion / extension settings
- ActiveSeam stitch

INDICATIONS

- Anterior Knee Pain
- Maltracking Patella
- Patellofemoral Osteoarthritis
- Chondromalacia

Product

Q Lok ROM - Closed Patella

Q Lok ROM - Closed Patella - Front Closure

Product

6166*

5166*

Note: To indicate size, replace * with corresponding number from the sizing chart on the back cover.



Adjustable Hinge for Q Lok ROM & CrossFire ROM

CrossFire ROM™ Knee Ligament Protection Featuring the new SlimLok™ Adjustable Hinge

The CrossFire ROM is designed to protect knees with ligament tears or sprains. The rigid uprights are slightly flared to contour to the leg's natural shape and provide stability so patients can move and be active with confidence. The thin hinge is easy to adjust for flexion and extension control during patient rehabilitation.

FEATURES

- Length: 13"
- Fits left or right
- Easily adjustable flexion / extension settings
- ActiveSeam stitch

INDICATIONS

- Ligament Sprains / Tears
- Meniscus Tears
- Knee Instability

Product

CrossFire ROM - Closed Patella

CrossFire ROM - Closed Patella - Front Closure

Product

6180*

5180*

Note: To indicate size, replace * with corresponding number from the sizing chart on the back cover.

Hinged Knee Skin™ Lightest hinged knee brace

Our new hinge is 25% lighter and twice as strong. Much thinner too, with sleek uprights (only 1.5mm) and smooth-motion hinge gears.

FEATURES

- *Offset hinge reduces pressure at the knee joint*
- *Strapping system minimizes migration*
- *Open patella*
- *Pull-on **or** front-closure design*
- *Fits left or right*

INDICATIONS

- *LCL Tear or Sprain*
- *MCL Tear or Sprain*
- *Meniscus Tear*

Product

Product	Product #
Hinged Knee Skin - Open Patella - Ultima™ Material	4130*
Hinged Knee Skin - Open Patella - Stratus™ Material	6130*
Hinged Knee Skin - Front Closure	5133*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Pull-on

Front-closure

Q Lok™ Dynamic Patella Traction System

The Q Lok has been shown to affect patellar positioning in weight-bearing MRI studies. Our unique design decreases pain by increasing surface contact area in the joint.

FEATURES

- *Addresses glide and tilt*
- *Closed patella*
- *Lightweight: 7 ounces*
- *Pull-on **or** front-closure design*
- *Fits left or right*

INDICATIONS

- *Patellar Tilt and Glide*
- *Patellar Subluxation*
- *Post Surgical Lateral Release*
- *Patellofemoral Osteoarthritis*

Product

Product	Product #
Q Lok Dynamic Patella Traction Brace - Stratus Material	6164*
Q Lok Dynamic Patella Traction Brace - Front Closure	5164*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Pull-on

Front-closure





Conforma™ Hinge

The perfect solution for patients who need extra stability, but also want something that's light and comfortable.

Lightweight

Each hinge weighs only one ounce, making the total weight for an entire brace a mere 6.75 ounces.

Flexible

Nylon plastic uprights conform perfectly to the shape of any leg without the need for any alterations by the fitter.

Low-Profile

The Conforma hinge is less than 4/10 of an inch thick and the uprights are less than 1/10 of an inch thick.

Visco Knee Skin™ with Conforma™ Hinge

Silicone ring reduces swelling around the joint.

FEATURES

- Ideal after knee arthroscopy
- Finger loops enable easy pull-on application
- Closed patella
- Fits left or right

INDICATIONS

- Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis

Product	Product #
Visco Knee Skin with Conforma Hinge	6137*
- Closed Patella - Stratus™ Material	

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

Q Brace™ with Conforma Hinge

Applies dynamic, lateral pressure to the patella.

FEATURES

- Open patella
- Pull-on **or** front-closure design
- Fits left or right

INDICATIONS

- Anterior Knee Pain
- Patellar Subluxation
- Lateral Patella Tracking
- Medial/Superior Patella Tracking
- Medial Patella Tracking
- Medial/Inferior Patella Tracking

Product	Product #
Q Brace w/ Conforma Hinge	4136*
- Open Patella - Ultima™ Material	
Q Brace w/ Conforma Hinge - Front Closure	5138*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

Patella Stabilizer with Conforma Hinge

Stabilizes the patella and provides general support.

FEATURES

- Foam buttress around patella for stability
- Open patella
- Fits left or right

INDICATIONS

- Patella Stabilization
- Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis

Product	Product #
Patella Stabilizer w/ Conforma Hinge	4132*
- Open Patella - Ultima Material	
Patella Stabilizer w/ Conforma Hinge	6132*
- Open Patella - Stratus Material	

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

Premium J™ with Conforma Hinge

Simple solution for patellofemoral pain.

FEATURES

- Adjustable, lateral patellar tension straps
- Closed patella
- Specify left or right

INDICATIONS

- Anterior Knee Pain
- Lateral Tracking Patella
- Patellar Subluxation
- Patellofemoral Knee Pain

Product	Product #
Premium J with Conforma Hinge	6138*L
- Closed Patella - Stratus Material - Left	
Premium J with Conforma Hinge	6138*R
- Closed Patella - Stratus Material - Right	

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

Q Brace™

This unique, multi-functional, patella-tracking brace offers multi-directional tracking to the patella.

FEATURES

- Achieves greater dynamic force as the knee is flexed
- Open patella
- Pull-on **or** front-closure design
- Fits left or right
- Available with a hinge for increased stability (see pg. 8)

INDICATIONS

- Anterior Knee Pain
- Patellar Subluxation
- Lateral Patella Tracking
- Medial/Superior Patella Tracking
- Medial Patella Tracking
- Medial/Inferior Patella Tracking

Product	Product #
Q Brace - Open Patella - Ultima™ Material	4161*
Q Brace - Front Closure	5163*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Pull-on Front-closure

Visco Knee Skin™

Silicone ring reduces swelling and provides passive massage to control edema.

FEATURES

- Ideal after knee arthroscopy
- Pull-on finger loops
- Closed patella
- Fits left or right
- Available with a hinge for increased stability (see pg. 8)

INDICATIONS

- Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis

Product	Product #
Visco Knee Skin - Closed Patella - Ultima Material	4170*
Visco Knee Skin with Straps - Closed Patella - Ultima Material	4171*
Visco Knee Skin with Straps - Closed Patella - Stratus™ Material	6171*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Q Baby™

Low-profile compression below the patella.

FEATURES

- Silicone pad may be placed over patellar tendon or tibial tuberosity
- Fits left or right

INDICATIONS

- Patellar Tendonitis
- Osgood-Schlatter Disease
- Patellofemoral Knee Pain

Product	Product #
Q Baby Patellar Tendon Strap	5165*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Premium J™

Simple solution for patellofemoral pain.

FEATURES

- Adjustable, lateral patellar tension
- Closed patella
- Specify left or right
- Available with a hinge for increased stability (see pg. 8)

INDICATIONS

- Anterior Knee Pain
- Lateral Tracking Patella
- Patellar Subluxation
- Patellofemoral Knee Pain

Product	Product #
Premium J - Stratus™ Material - Left	6163*L
Premium J - Stratus™ Material - Right	6163*R

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Patella Stabilizer

Stabilizes the patella and provides general support.

FEATURES

- Foam buttress ring around patella for stability
- Open patella
- Fits left or right
- Available with a hinge for increased stability (see pg. 8)

INDICATIONS

- Patella Stabilization
- Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis
- Patellar Tendonitis

Product	Product #
Patella Stabilizer - Open Patella - Ultima™ Material	4125*
Patella Stabilizer - Open Patella - Stratus™ Material	6125*

Note: To indicate size, replace * with corresponding number from the sizing chart on the back cover.



No Straps

With Straps

Standard Knee Skin™

Compression for minor knee pain and instability.

FEATURES

- Available in open or closed patella
- Fits left or right
- Available with straps

INDICATIONS

- Mild Knee Sprains and Strains
- Knee Joint Swelling
- Meniscus Tear
- Osteoarthritis

Product	Product #
Standard Knee Skin - Open Patella - Ultima Material	4110*O
Standard Knee Skin - Closed Patella - Ultima Material	4110*C
Standard Knee Skin with Straps - Open Patella - Ultima Material	4111*O
Standard Knee Skin with Straps - Closed Patella - Ultima Material	4111*C

Note: To indicate size, replace * with corresponding number from the sizing chart on the back cover.

Under Skin™

Provides compression and prevents skin irritations while wearing a rigid brace.

FEATURES

- Thin, breathable compressive material
- 20" in length
- Decreases muscle oscillation
- Limits migration of a rigid brace
- Latex-free and hypo-allergenic

INDICATIONS

- Protects Skin from Rigid Braces
- Post-Workout Recovery
- Muscle Injury
- Mild Knee Swelling

Product	Product #
Underskin	4140*

Note: To indicate size, replace * with corresponding number from the sizing chart on the back cover.

BIO SKIN Materials Are:

Compressive

Hypo-Allergenic

Neoprene-Free

Latex-Free

SkinLok™ Anti-Slip



Why choose *Baja*™?

Pain Relief – Comfort – Ease of Use – High Value

Pain Relief

- **Clinical validation**

A recent study using the bilateral cinching mechanism of the Baja showed “Patients’ average VAS pain score dropped by 41.2% immediately upon donning the device...”

*“Prospective Study of Lumbar Back Brace In an Interventional Pain Practice.” Rosenthal, Richard M.D. & Spencer, Shawn, BS. Practical Pain Management, (July/August 2011)

- **Vectored compression**

Delivers targeted pain relief directly to the lower spine

- **Integrated hot/cold gel pads**

Easily slides into the brace pouch for additional pain relief and comfort



Hot/Cold Gel Pads



Breathable Fabric

Comfort

- **Breathable spacer fabric**

Allows heat dissipation and quick evaporation

- **Contoured waist belt**

Fits comfortably around the rib cage and iliac crest

- **Lightweight**

Weighs less than 2 pounds

- **Slim profile**


Sleek design conforms to the body



Ergonomic Hand Pulls & Contoured Fit

Ease of Use

- **Bilateral cinching**
Bilateral cinching mechanism doubles the compression provided with equal force
- **Mechanically-advantaged rollers**
Vertical rollers reduce the force needed to apply compression
- **Ergonomic hand pulls**
Substantial handles are easy to pull, even for arthritic patients
- **Easy-release Velcro**
Advanced Velcro provides firm grip and easy release



Bilateral Cinching

High Value

- **Available in pre-sized or universal, trim-to-fit model**
Choose the convenience of dispensing off-the-shelf or minimizing inventory with trim-to-fit model
- **High patient acceptance rate**
High acceptance rate by patients maximizes clinic revenue
- **1-year warranty**
Cropper Medical warranties all Bajas against defects for 12 months from purchase date
- **Custom logos**
Increase your clinic's exposure by adding your logo (call for details)




Mechanically-Advantaged Rollers



Custom Logo


How the Baja Back Brace Line is Different

Baja – Posterior Compression



The Baja has a bilateral lacing system with two banks of rollers positioned on the sides of the torso. When the patient engages the cinching mechanism, the posterior panel is drawn forward, focusing the compression on the soft tissues that support the spine.

Other Brands – Lateral Compression



Other cinch-style back braces have a *single* posterior lacing system with the laces positioned *across* the spine. When a patient engages this style of brace, the two sides pull together and the compression is focused on the *sides* of the torso rather than the spine.



Baja™ 627

FEATURES

- Posterior panel provides support from **L-5 to L-1**
- Two hot/cold gel pads
- Easy-to-pull bilateral cinching system

INDICATIONS

- Herniated Discs
- Degenerative Disc Disease
- Lumbar Sprains and Strains
- Spondylolisthesis
- Osteoarthritis

Product	Product #
Baja 627	3750*
Baja 627 - Trim to Fit - Universal Size	37504U

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Baja™ 631

FEATURES

- Posterior panel provides support from **sacrum to T-9**
- Two hot/cold gel pads
- Easy-to-pull bilateral cinching system

INDICATIONS

- Herniated Discs
- Degenerative Disc Disease
- Spondylolisthesis
- Osteoarthritis
- Radicular Pain

Product	Product #
Baja 631	3770*
Baja 631 - Trim to Fit - Universal Size	37704U

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Baja™ 637

FEATURES

- Posterior panel provides support from **sacrum to T-9**
- Two hot/cold gel pads
- Easy-to-pull, bilateral cinching system
- Removable side panels to limit coronal motion
- Sagittal control

INDICATIONS

- Herniated Discs
- Degenerative Disc Disease
- Spondylolisthesis
- Osteoarthritis
- Radicular Pain
- Lumbar Spinal Stenosis
- Post-Surgical Stabilization

Product	Product #
Baja 637	3760*
Baja 637 - Trim to Fit - Universal Size	37604U

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Back Flex™

A simple and effective treatment for a variety of causes of lumbar pain.

L0627
L0642

HCPs

FEATURES

- Two hot/cold gel pads
- Cinching bilateral pulley system

INDICATIONS

- Herniated Discs
- Degenerative Disc Disease
- Lumbar Sprains & Strains
- Spondylolisthesis
- Osteoarthritis

Product

Back Flex

Product

3730*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Back Plus™

The cinching mechanism allows a patient to easily apply circumferential compression.

L0631
L0648

HCPs

FEATURES

- Two hot/cold gel pads
- Cinching bilateral pulley system
- Shoulder straps

INDICATIONS

- Compression Fractures
- Post-Vertebroplasty
- Post-Kyphoplasty
- Post-SCS Implant
- Post-IT Pump Implant
- Post-Fusion

Product

Back Plus

Product

3740*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Back Skin™

Medium-level compression for muscle strains and minor back pain.

FEATURES

- Narrow cut in front
- Removable double-cinch belt
- Optional lumbar accessories

INDICATIONS

- Lumbar Strains
- Lumbosacral Strains
- Lumbago
- Lower Back Pain

Optional Lumbar Accessories



Lumbar Pad



Flexible Support

Product

Standard Back Skin

Standard Back Skin with Flexible Support

Standard Back Skin with Lumbar Foam Pad

Standard Back Skin with Flexible Support & Lumbar Pad

Product

4710*

4734*

4720*

4724*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

TriLok®

External Ligament System™

The only brace to offer Medial & Lateral Unloading

VERSATILITY:

Can be applied for medial or lateral pathologies and worn on either the left or right foot.

CONTROL:

Two independent straps stabilize the midfoot and hindfoot.

COMFORT:

Low profile allows for normal ankle flexion and extension. Wear barefoot, in shoes or sandals.

INDICATIONS

- Posterior Tibial Tendon Dysfunction
- Plantar Fasciitis
- Peroneal Tendinopathy
- Ankle Sprains
- Subtalar & Ankle Arthritis

HCPCS

L1902

Product

TriLok

Product

5360*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



TriLok® Unequaled Versatility



The TriLok has unparalleled versatility, delivering targeted unloading to medial and lateral structures. Positioned medially, the FootLok™ Strap supports the plantar fascia or the posterior tibial tendon. Positioned laterally, the FootLok Strap supports and protects the lateral ligaments and the peroneal tendons. The TriLok is effective with almost any type of footwear or with no footwear at all.

Medial Unloading:

Plantar Fasciitis – Posterior Tibial Tendon Dysfunction – Deltoid Sprains – Subtalar Malalignment



Lateral Unloading:

Lateral Sprains – Peroneal Tendon Damage – Ankle Arthritis



The AFTR™ Line combines stabilization with high-level compression unique to Bio Skin®. Following trauma, AFTR braces help:

- Reduce Pain
- Improve Proprioception
- Increase Range of Motion
- Reduce Recovery Time

INDICATIONS

- Peroneal / Post-Tib Tendon Repair
- Lateral Ankle Repair
- Fusions: Midfoot, Hindfoot, Ankle
- Total Ankle Replacement
- Ankle Sprains
- Bunionectomy
- Forefoot Procedures
- Hindfoot Arthritis

AFTR DC with Gel™ (Dorsal Closure)

FEATURES

- Accommodates changes in foot volume
- Easy to apply during early, post-op phase
- Includes rigid stirrup strap
- Includes medial and lateral silicone gel pads
- Trimmable in length
- Available with compression wrap

Product	Product #
AFTR DC	5386*
AFTR DC with Compression Wrap	5386*-CW

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

AFTR with Gel™ (Pull-on)

FEATURES

- Pull-on design
- Includes medial and lateral silicone gel pads behind the malleoli
- Includes rigid stirrup strap
- Trimmable in length
- Available with compression wrap

Product	Product #
AFTR with Gel	5383*
AFTR with Gel with Compression Wrap	5383*-CW

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

AFTR™ (Pull-on)

FEATURES

- Pull-on design
- Trimmable in length
- Includes rigid stirrup strap
- Available with compression wrap

Product	Product #
AFTR	5380*
AFTR with Compression Wrap	5380*-CW

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Compatible w/
AFTR braces.



Shown with [Calibrated Midfoot Compression Wrap](#) (sold separately).
Also compatible with [AFTR DC with Gel](#), [AFTR with Gel](#) and [AFTR](#).

HCS™ (Hallux Control Strap™)

29550

CPT

Splinting - Maintains the hallux in a neutral position during critical soft-tissue healing after bunionectomy.

Stretching - Stretches extensor tendons after primary healing to improve ROM and allow normal gait.

Strengthening - Provides resistance for patients to strengthen FHL for strong toe-off.

Note: Must be worn with any AFTR brace, or the Calibrated Midfoot Compression Wrap.

Product	Product #
HCS 10-Pack	9534*†

Compatible w/
AFTR braces.



Shown with [Calibrated Midfoot Compression Wrap](#) (sold separately).
Also compatible with [AFTR DC with Gel](#), [AFTR with Gel](#) and [AFTR](#).

Weil Osteotomy Strap™

29550

CPT

Prevents dorsal drift (floating toe) after Weil Osteotomy by holding the toe in a neutral or plantar-flexed position.

May be worn inside a walking boot, post-op shoe, or normal footwear.

May be used non-operatively to treat metatarsalgia or plantar plate tears.

Note: Must be worn with any AFTR brace, or the Calibrated Midfoot Compression Wrap.

Product	Product #
Weil Osteotomy Strap 10-Pack	95372



Compatible with the [Hallux Control Strap](#)
and [Weil Osteotomy Strap](#).

Calibrated Midfoot Compression Wrap™

29540

CPT

Provides an attachment point for the Weil Digit Strapping System without covering the ankle.

Made of hypo-allergenic, non-slip material to prevent migration.

Easily fits inside normal footwear.

Product	Product #
CMC Wrap 10-Pack	95362



DP2 - 6 1/2" DP3 - 8 1/2"

DP2™ & DP3™ Cock-Up Wrist

L3908

HCPs

FEATURES

- Low profile and comfortable
- Soft, fleece material against skin
- Allows full use of fingers and thumbs
- Removable/bendable dorsal and palmar stays
- Specify left or right

INDICATIONS

- Carpal Tunnel Syndrome
- Arthritis
- Wrist Sprain
- Bursitis

Product	Product #
DP2 Cock-Up Wrist Splint 6 1/2" - Left	5423*L
DP2 Cock-Up Wrist Splint 6 1/2" - Right	5423*R
DP3 Cock-Up Wrist Splint 8 1/2" - Left	5424*L
DP3 Cock-Up Wrist Splint 8 1/2" - Right	5424*R

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Wrist/Thumb Spica

L3807

HCPs

FEATURES

- Controls motion of the CMC and MP joints
- Cock-up wrist support
- Soft, fleece material against skin
- Specify left or right

INDICATIONS

- De Quervains
- Scaphoid Fracture
- Bursitis
- Sprains and Strains

Product	Product #
Wrist/Thumb Spica - Left	5452*L
Wrist/Thumb Spica - Right	5452*R

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Thumb Spica

L3923

HCPs

FEATURES

- Controls motion of the CMC and MP joints
- Low profile
- Soft, fleece material against skin
- Fits left or right

INDICATIONS

- De Quervains
- Game Keeper's Thumb
- Bursitis
- Sprains and Strains

Product	Product #
Thumb Spica	5450*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Boomerang

A4466

HCPs

FEATURES

- Universal size
- Easy attachment
- Soft, fleece material against skin
- Fits left or right

INDICATIONS

- General Wrist Support
- Carpal Tunnel Syndrome
- Arthritis
- Wrist Sprain
- Wrist Pain
- Bursitis

Product	Product #
Boomerang Wrist Wrap	54200



DP2 Cock-Up Wrist



Standard Elbow Skin™

FEATURES

- Full elbow coverage
- Comfortable, low-profile fit

INDICATIONS

- Bursitis
- Joint Swelling
- Biceps/Triceps Strains

Product	Product #
Standard Elbow Skin - Ultima Material	4210*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Tennis Elbow Skin™

FEATURES

- Silicone pad fits over painful area
- Full elbow coverage

INDICATIONS

- Medial Epicondylitis (Golfer's Elbow)
- Lateral Epicondylitis (Tennis Elbow)

Product	Product #
Tennis Elbow Skin - Ultima Material	4220*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Tennis Elbow Band™

FEATURES

- Silicone pad fits over painful area
- Fits comfortably below elbow

INDICATIONS

- Medial Epicondylitis (Golfer's Elbow)
- Lateral Epicondylitis (Tennis Elbow)

Product	Product #
Tennis Elbow Band - Ultima Material	4250*
Tennis Elbow Band - 2S Material	5250*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

Standard Elbow Skin





Calf Skin™

FEATURES

- Controls muscle movement
- Reduces muscle vibration

INDICATIONS

- Contusion
- Shinsplints
- Calf/Shin Pain

Product

Calf Skin - Ultima Material

Product

4630*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Optional Cinch Strap



Thigh Skin™

FEATURES

- Gives optimal compression to quads and hamstrings
- Optional Cinch Strap for quad and hamstring injuries

INDICATIONS

- Quad Contusion
- Thigh Strain

Product

Thigh Skin - Ultima Material

Thigh Skin with Cinch Strap - Ultima Material

Product

4610*

4611*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.



Optional Groin Wrap



Shorts in white for contrast only.

Compression Shorts

FEATURES

- Highly Compressive
- Muscle movement control to improve recovery time
- Increases proprioception
- Optional Cinch Strap and Groin Wrap for isolated compression

INDICATIONS

- Quad Contusion
- Hip Flexor Strain
- Ischiocapsular/Iliofemoral/ Thigh Strain

Product

Compression Shorts - Ultima Material

Compression Shorts with Groin Wrap - Ultima Material

Product

4620*

4621*

Note: To indicate size, **replace** * with corresponding number from the sizing chart on the back cover.

Calf Skins




► PREMIUM MATERIALS

We've designed premium materials that are hypo-allergenic, neoprene-free, latex-free, breathable, light, thin, and strong. These special materials form the foundation upon which we build our entire product line and allow us to provide your patients with the best and most comfortable braces on the market.

Ultima™

Ultra-Thin, High Compression



Ultima is our original, tri-laminate material. Ultima is amazingly thin and provides high-level compression with powerful Lycra and a polyurethane film. The film allows moisture vapor transfer so perspiration evaporates, facilitating the body's natural temperature regulation system. The Ultima material is used widely in our product line.

Stratus™

Extra Breathability, Forgiving Fit



Stratus is slightly thicker than Ultima for people who like the feeling of extra support. Stratus is a laminated material with high-power Lycra on the outsides and a porous membrane in the middle. The more elastic membrane in Stratus allows it to provide more dynamic compression and a more forgiving fit. The Stratus is used primarily in some of our knee products.

BIO SKIN Materials Are:

Compressive

Hypo-Allergenic


Neoprene-Free

Latex-Free

SkinLok™ Anti-Slip

2SL™


Front Closure Technology



2SL features a Velcro-compatible surface to provide compression and flexibility in anchoring and closures. The 2SL material is used primarily in our Front-Closure Knee Skins to provide compression and allow total adjustability.

2S™

Soft, Micro-Fleece Lining




2S material is similar to our Ultima material, but the skin-side of the fabric is made of moisture-wicking micro-fleece. The fleece is extra soft against the skin and is treated with an anti-microbial finish to prevent odor. The 2S material is used in several of our wrist and ankle products.

NEW

Cirrus™


Single Layer, High Performance Stretch



Cirrus is made of a single layer of high-performance Lycra. Cirrus is incredibly thin and perfectly breathable, but still provides a moderate level of compression. The Cirrus material is used in certain knee braces to prevent migration and limit bunching behind the knee.

5200™

Compressive Closure Technology



5200 provides moderate compression with soft, brushed Nylon/Lycra. It's soft against the skin and is Velcro-compatible to allow attachment of straps and accessories on the AFTR line of ankle braces.

Knee

Knee Products					
Circumference with knee fully extended					
Size	*PN	Inches		Centimeters	
		mid – patella	7" above	mid – patella	18cm above
XS	0	12 – 13	16 – 18	31 – 33	41 – 46
S	1	13 – 14	18 – 20	33 – 36	46 – 51
M	2	14 – 16	20 – 22	36 – 41	51 – 56
L	3	16 – 17	22 – 25	41 – 43	56 – 64
XL	4	17 – 19	25 – 28	43 – 46	64 – 71
2XL	5	19 – 21	28 – 32	46 – 54	71 – 81

Thigh

Mid – thigh circumference				
Size	*PN	Inches	Centimeters	
XS	0	16 – 18	41 – 46	
S	1	18 – 20	46 – 51	
M	2	20 – 22	51 – 56	
L	3	22 – 25	56 – 64	
XL	4	25 – 28	64 – 71	
2XL	5	28 – 32	71 – 81	

Spine & Shorts

Circumference between hip and top of pelvis.			
Size	*PN	Inches	Centimeters
Baja™			
S	1	26 – 30	66 – 76
M	2	30 – 36	76 – 91
L	3	36 – 44	91 – 112
XL	4	44 – 54	112 – 137
UNIV	4U	26 – 54	66 – 137
Back Flex & Back Plus			
1	1	28 – 34	71 – 86
3	3	34 – 44	86 – 112
5	5	44 – 56	112 – 142
Back Skin			
XS	0	24 – 26	61 – 66
S	1	26 – 30	66 – 76
M	2	30 – 34	76 – 86
L	3	34 – 38	86 – 97
XL	4	38 – 42	97 – 107
2XL	5	42 – 46	107 – 117
Compression Shorts			
S	1	28 – 32	71 – 81
M	2	32 – 36	81 – 92
L	3	36 – 40	92 – 102
XL	4	40 – 44	102 – 112
2XL	5	44 – 48	112 – 122

Elbow

Circumference of largest part of forearm				
Size	*PN	Inches	Centimeters	
XS	0	8 – 9	20 – 23	
S	1	9 – 10	23 – 26	
M	2	10 – 11	26 – 28	
L	3	11 – 12	28 – 31	
XL	4	12 – 13	31 – 33	
2XL	5	13 – 14	33 – 36	

Calf

Largest part of calf circumference				
Size	*PN	Inches	Centimeters	
S	1	11 – 12	28 – 30	
M	2	12 – 14	30 – 35	
L	3	14 – 16	35 – 41	
XL	4	16 – 18	41 – 46	
2XL	5	18 – 20	46 – 51	

Call for custom products and extended sizing.

Continued on Next Page

Indicate Sizing When Ordering

Product#	Size#	To indicate size, replace * with corresponding *PN (Part Number) from the sizing chart. Custom products and extended sizing is available. For more information, please call: 1-800 -265-7126
XXXXX	*	



ON-LINE ORDERING: www.DunbarMedical.com

SIZING CHART

Foot & Ankle

By shoe size					
TriLok™					
Size	*PN	Men's	Women's	Europe	UK
S	1	5 – 6.5	6 – 7.5	38 – 39	4 – 5.5
M	2	7 – 9.5	8 – 11	40 – 42	6 – 8.5
L	3	10 – 12.5		43 – 45	9 – 11.5
XL	4	13 – 16		46 – 49	12 – 16
HCS - Hallux Control Strap					
S	1	5 – 10	6 – 11	38 – 43	4 – 9
L	3	10.5 – 16	11.5+	44 – 50	9.5 – 16

Circumference from heel to top of the foot – back to heel			
Size	*PN	Inches	Centimeters
AFTR Products			
XS - S	1	11 – 13	27.9 – 33
M - L	3	12.5 – 14	31.8 – 35.6
XL - 2XL	5	13.5 – 15	34.3 – 38.1
Standard & Visco			
XS	0	10 – 11	25 – 28
S	1	11 – 12	28 – 31
M	2	12 – 13	31 – 33
L	3	13 – 14	33 – 36
XL	4	14 – 15	36 – 38
2XL	5	15 – 16	38 – 41

Wrist

Circumference of extended wrist			
Size	*PN	Inches	Centimeters
XS - S	1	4 – 6	10 – 15
M - L	3	6 – 8	15 – 20
XL - 2XL	5	8 – 10	20 – 25

Thumb Spica			
Length of thumb			
Size	*PN	Inches	Centimeters
XS - M	2	3.25 – 4.25	8 – 11
L - 2XL	5	4.25 – 5.25	11 – 14

Four Easy Ways to Order



ON-LINE

NEW!

www.DunbarMedical.com

Easy, Fast, Secure, 24/7

- Place Orders
- Track Shipments
- Application Videos
- Product Reviews



EMAIL

Sales@DunbarMedical.com

Email orders 24/7

Use Subject Line: "Bio Skin Order"



PHONE

1-800- 265-7126

Our Customer Care team is ready to help. Monday - Friday, 9am - 5pm Eastern time.



FAX

1-800-263-5607 (Fax orders 24/7)

Call for custom products and extended sizing.

Continued on Previous Page

Indicate Sizing When Ordering

Product#	Size#	
XXXXX	*	To indicate size, replace * with corresponding *PN (Part Number) from the sizing chart. Custom products and extended sizing is available. For more information, please call: 1-800- 265-7126



BioSkin™ 5285 Tomken Road, Unit 1, Mississauga, ON L4W 3N8



t: 800-265-7126

f: 800-263-5607

www.DunbarMedical.com



Working Together...To Help Others



Terms & Conditions

NET 30 ACCOUNTS

1. All "net 30" invoices are payable within 30 days of shipping date. A charge of 2% per month will be added to all invoices outstanding for more than 30 days.
2. A handling charge of \$25 will be assessed for all "NSF" cheques.
3. All new accounts will be shipped on a prepaid or COD basis until credit approval is obtained.

PAYMENTS

1. Pay oldest invoice first.
2. Include invoice or statement stub with payment.
3. Discount prices are available only to customers with no past due accounts.

SHIPPING POLICY

1. Delivery to carrier constitutes delivery to purchaser. Claims for missing merchandise must be filed by the recipient within three (3) business days of the expected date of delivery.
2. Orders of \$500 or more before taxes within Ontario (\$800 or more outside Ontario) are exempt from ground shipping charges. Orders that include liquids, gels, Renfrew Professional Trainers Tape, Hot/Cold Packs and ice bags are not eligible for this offer.

RETURN AND REPLACEMENT POLICY

1. Merchandise may not be returned without our approval. Please obtain a return authorization number before returning an item.
2. All returns are subject to inspection.
3. A 25% restocking fee will be applied to all returns, unless the item returned is being exchanged for another of the same type but a different size or color.
4. Shipping for all authorized returns and for exchanged items to be paid by the customer.
5. Include a note explaining the reason for return, purchase date and invoice number.
6. Products purchased with initial stocking order may be exchanged within three (3) months without a restocking fee.
7. Except for initial orders as stated above, products must be returned within 30 days.
8. Repairs must be clean and dry or they will be returned to you without repair.
9. Custom or special sized products may be returned for modification only, not for credit or refund.
10. ALL RETURNS MUST BE IN SALEABLE CONDITION AND IN ORIGINAL PACKAGING.

PRICE CHANGES

1. **All prices are subject to change without notice.**

