

## FEATURES

- *Low-Profile and Ultra Lightweight*
- *Rounded Sole for Improved Maneuverability*
- *Pre-Shaped Ergonomic Frame*
- *Automatically Adjusts to Gait Patterns*
- *New and improved pump system*
- *Absorbs 50% of heel strike shock*
- *HCPCS: L4361*



# PNEUMATIC WALKING BOOT

Low Profile — Ultra Lightweight Frame



By shoe size			
Size	PN	Men's	Women's
XS	33100	2 - 4	3 - 6
S	33101	4 - 7	6 - 8
M	33102	7.5 - 10.5	8 - 11
L	33103	10.5 - 12.5	11 - 113
XL	33104	12.5+	13+





# BioSkin™ Pneumatic Walking Boot

A breakthrough low-profile design providing higher maneuverability, stability and heel relief.



## Pain Relief

- Flared contoured struts to fit a wide range of patients
- Comprehensive leg and ankle protection
- Post-surgical Immobilization of the foot and ankle
- Decrease pain and protect foot for enhanced rehabilitation

## Comfort

- Pre-relieved heel contact area to increase comfort and safety
- Low-profile, rounded sole for enhanced maneuverability
- ShockPod™ system absorbs over 50% of force from heel strike
- Wide foot bed for improved comfort and stability

## Ease of Use

- Adjustable pneumatic pump with easy “push-button” release
- Heel height of a traditional shoe for ease of ambulation
- Fully adjustable strapping for sizing accommodations

## Indications

- Ankle Fractures
- Post-Operative Immobilization
- Post-Operative Rehabilitation
- Grade 2 and 3 Ankle Sprains
- Soft Tissue Injuries



Apply the boot liner



Place foot into boot frame



Secure straps



Inflate pneumatic pump



240 East Hersey St., Suite 2, Ashland, OR 97520-5202

t: 800-541-2455 f: 541-482-2119 Pro.BioSkin.com

©2016 Cropper Medical, Inc.  
4/2016 PN: 12509A